

Agilent 2100 Bioanalyzer – Sample Submission Form

Client Details

Name: _____	
Phone: _____	
Email: _____	
Billing Address: _____ _____ _____	Purchase Order Number: _____
	Authorisation (please sign): _____

Sample Details

NB. Samples submitted for analysis using the RNA Nano chip should be in a volume of 3 uL at 50-500 ng/uL.

RNA source: Eukaryote Prokaryote

Sample	Name	Concentration (ng/uL)	Volume (uL)
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			
11			
12			